


Manage The Entire Endoscope Reprocessing Loop

ScopeTrac® Advanced manages every endoscope through the entire procedural cycle – from bedside to storage.

- Track endoscopes electronically
- Process up to four scopes simultaneously with safety stops and scope messages
- Paperless reprocessing documentation
- Guided workflows and IFU access
- Comprehensive data and reporting capabilities

- Interfacing capabilities with OR schedulers and reprocessing machines
- Enhance workflow efficiency with hands-free optimization
- Scope sequence warnings
- Request maintenance while processing scopes



Decontamination

ScopeTrac Advanced identifies endoscopes through barcode scanning technology, or by letting the user manually select the endoscopes.

Clean Room

ScopeTrac Advanced displays scope-specific customer-generated and/or manufacturer instructions related to high-level disinfection (HLD) and sterilization for the user to follow step by step.

Clean Storage

As endoscopes and their accessories are scanned into storage, assets currently available can be accessed from any workstation by authorized users.

Procedure

By scanning the endoscope barcode on the asset itself, ScopeTrac Advanced allows the user to link a specific scope to a specific patient case and store this information as part of the scope history.



ScopeTrac Advanced Reporting

ScopeTrac Advanced provides comprehensive, real-time management reports that can be run on-demand and provide visibility at both summary and detailed levels. Endoscope reprocessing report customization options include the ability to:

- Specify any combination and sequence of data fields
- Filter data based on user-selected values
- Hide data fields not required in the customer report
- Set up exception reporting that can alert users via email of issues requiring attention
- Extrapolate reprocessing histories to include testing, cleaning, point-of-use treatment with user and real-time data

